

A person wearing a white bucket hat, a light-colored jacket, blue jeans, and yellow boots is fishing on a wooden dock. They are holding a fishing rod with a reel. A portable oxygen concentrator (PHILIPS SimplyGo) is attached to their waist with a shoulder strap. The device is olive green and tan. In the background, there is a body of water and trees. A tackle box is visible on the dock to the left.

**PHILIPS**  
RESPIRONICS

*SimplyGo*

Portable oxygen  
concentrator

Go with  
**SimplyGo**

A great choice in  
continuous flow portable  
oxygen concentrators  
weighing ten pounds  
or less



# Good to go with SimplyGo

Whether it's a walk across the street, a ride across town, or a trip across the country, Philips Respironics has made it easier than ever for nearly all oxygen users to get around with the SimplyGo portable oxygen concentrator (POC).

SimplyGo is truly a new breed of lightweight and attractive portable oxygen concentrator with continuous flow capability.

At just ten pounds, SimplyGo can deliver up to four times more oxygen than some other POCs.



## Fits your lifestyle

SimplyGo's non-clinical style and ease-of-use make it appealing to oxygen users. An attractive and highly functional carrying case allows users to either carry SimplyGo over the shoulder or place it on the easy-to-pull, easy-to-store cart.

## Redefines ruggedness

SimplyGo was built to set new standards of ruggedness in order to withstand the bumps and bruises of daily activities and deliver reliable performance. A sturdy cart with a locking adjustable handle and oversized wheels make SimplyGo ready for many of life's obstacles.

## Powered for convenience

One of the keys to SimplyGo's portability is its lightweight and compact lithium ion battery that is immediately accessible and can be easily changed. Standard AC and DC power cords add convenience and flexibility.

---

The SimplyGo system comes with one lithium ion battery. Extra batteries fit easily into the zippered pouch or handy accessory bag. Want more time? Simply take along more batteries.





## Is SimplyGo right for me?

The SimplyGo portable oxygen concentrator is intended for prescription use by patients requiring high concentrations of oxygen on a supplemental basis. It may be considered for continuous use in home, institutional, and travel/mobile environments. SimplyGo is not intended to be life-sustaining or life-supporting. And is not intended for pediatric use. SimplyGo may not be for all oxygen therapy needs. You should talk to your doctor to determine if it is right for you.



Because SimplyGo has more than twice the oxygen output of some other POCs of its size, it is able to keep up with users' portable oxygen needs today and tomorrow. SimplyGo offers both continuous flow (up to 2 liters per minute) and pulse-dose delivery (settings 1 to 6) while weighing only 10 pounds.

---

Visit [www.philips.com/simplygo](http://www.philips.com/simplygo)  
or ask your homecare provider for details.



© 2017 Koninklijke Philips N.V. All rights are reserved. Trademarks are the property of Koninklijke Philips N.V. or their respective owners. Respironics Inc. is a Philips company.

Specifications are subject to change without notice.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Caution: US federal law restricts these devices to sale by or on the order of a physician.

1010 Murry Ridge Lane, Murrysville, PA 15668  
800-345-6443 • 724-387-4000

[www.philips.com/respironics](http://www.philips.com/respironics)

KR 03/20/17 MCI 4107902 PN 1134358 R00